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Setting The Action Zones

The zone markers may be set with a stylus by the healthcare professional to give the personal boundaries for implementing the asthma action plan.

How To Use The Peak Flow Meter

- Slide the pointer to the bottom of the scale
- Hold the device horizontally in front of your mouth
- Inhale as deeply as possible
- Seal your lips firmly around the mouthpiece
- Blow out for a second or longer as HARD and as FAST as possible
- Take care that the tongue or teeth do not obstruct the airflow. A 'spitting' action or coughing will give false readings
- Take note of the reading and slide the pointer back to the bottom of the scale
- Repeat this three times and record the highest reading
- Take action according to your personal asthma management plan

Cleaning & Disinfecting

Use of an ordinary alcohol wipe (IPA 70-90%) is recommended with special attention to the mouthpiece area. Visually contaminated devices should be cleaned and disinfected prior to use.

Single Person Home Use

Clean at least weekly. Replace after 3 years.

Clinic Use (requires mouthpiece adaptor)

Between each subject: Clean, disinfect, and fit a new disposable Safeway or BVF mouthpiece. It is recommended that peak flow meters used in clinic be replaced every year or after 650 patients whichever is earlier.

Warnings & Advisory Notices

Note: Please read all this information before using this Peak Flow Meter.

- This Vitalograph device is intended to help you assess how well your asthma is under control, by measuring your peak flow which can then be monitored over time
- Never attempt to dismantle the unit. This can cause faulty Peak Expiratory Flow scores
- If used at home symptoms must take precedence over device measurements
- If you think that the device is not reading correctly, advise your healthcare professional immediately
- If the device is used for longer than it's specified life, the accuracy of the device may deteriorate
- A damaged Peak Flow Meter should always be replaced
- Store in a clean dry place

Warranty

Your Vitalograph device is guaranteed for one year*.

* Excepting accidental / transit damage or inappropriate use of the device

CE Notice

Marking by the symbol CE indicates compliance of this device to the Medical Devices Directive of the European Community. Such marking is required that the Vitalograph device meets or exceeds the referenced technical standards.

FDA Notice

Caution: Federal Law restricts this device to sale by, or on the order of a physician.

Declaration Of Conformity

Product: Vitalograph Model 4300
Vitalograph hereby ensures and declares that the above product associated with this user manual, is designed and manufactured in accordance with the following QMS regulations and standards:

- European Medical Devices Directive 93/42/EEC, as amended. This device is classified as 1m per Annex IX of the MDD also meets the provisions of the Essential Requirements, Annex I, via compliance with Annex II of the Medical Devices Directive as per Article 11, section 3a, excluding point 4 of Annex II;
 - Canadian Medical Device Regulation (CMDR SOR/98-282);
 - FDA Quality System Regulation (QSR) 21 CFR 820;
 - EN ISO 13485:2003. Medical devices. Quality management systems. Requirements for regulatory purposes
- Certifying Body (For 93/42/EEC and CMDR): British Standards Institute, Certificate Nos. CE 00772, CE 85553, MD 82182, FM 83550.

Technical Specifications

Product: Peak Flow Meter

Model Number: 4300

Material: Recyclable ABS plastic

Accuracy: ± 10 L/min or $\pm 10\%$ of the reading

Repeatability: ± 10 L/min or $\pm 5\%$ of the reading

Highest resistance to flow:

0.00377kPa/L/min @ 720 L/min

Measurement Range: 50-800 L/min BTPS

Operating Temperature Range:

ATS/ERS Limits: 17-37°C. Design Limits: 10-40°C

Storage Conditions:

Temperature: 0-50°C. Relative Humidity: 10%-95%

Peak Flow Meter performance standards:

EN ISO 23747:2009; ATS/ERS:2005; AS/NZS4237;

Frequency: Profile A/B difference less than 15 L/min/12%

(Annex B; EN ISO 23747:2009)

Medical Safety Standard: Medical Devices Directive

93/42/EEC1993 as amended

Designed & Manufactured Under: ISO 13485:2003, FDA

21CFR820, CMDR & JPAL

Signed on behalf of Vitalograph (Ireland) Ltd.

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B. R. Garbe, Direktor

Vitalograph
Your respiratory partner

Peak Flow Meter

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
the difference is clear


CE Peak Flow Meter
0086 Ref. No. 43985
EN ISO 23747
Model 4300 S 10991691 439856

Manufacturer: Vitalograph, Ennis, Ireland

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

Panel 1 shows a hand pointing to a PEF device. To the left is a circular icon with three stylized people and a plus sign. To the right are four status icons: a green smiley face (top right), a yellow neutral face (middle right), a red sad face (bottom right), and a vertical bar indicator (left) with yellow, red, and green segments.


Panel 2 is labeled with a '1' in a circle. It shows a hand holding the PEF device with the mouthpiece open. A red arrow points to the right. An inset shows the internal filter mechanism.

Panel 3 is labeled with a '2' in a circle. It shows a hand holding the PEF device to a person's mouth. A curved arrow indicates the direction of the inhaled air. An inset shows the hand holding the device.

Panel 4 is labeled with a '3' in a circle. It shows a hand holding the PEF device away from a person's mouth. A curved arrow indicates the direction of the exhaled air. Below the hand, the text 'x3' is written. An inset shows the hand holding the device.

Panel 5 contains recycling instructions. The top row shows a circular icon with three people and a plus sign, followed by a PEF device, a plus sign, a disassembled device, a number '2' inside a circle, and a recycling symbol. The bottom row shows a disassembled device, a plus sign, another disassembled device, an equals sign, and a fully assembled device.

PEF (l/Min)	800	750	700	650	600	550	500	450	400	350	300	250	200	150	100	50

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Gezeichnet im Namen von Vitalograph (Ireland) Ltd.